

Aesculap Scope Container

for the STERRAD® sterilisation system



Aesculap Scope Container for the STERRAD® sterilisation system

The Aesculap Scope Container combined with the STERRAD® sterilisation systems

STERRAD®100S



STERRAD®200

STERRAD®NX

STERRAD®100NX




is the perfect combination for protecting your valuable optical equipment.



As regards the compatibility of the items to be sterilised with the STERRAD® procedure, please observe the instructions for use of the respective manufacturer.

The following Aesculap products are compatible with the STERRAD® sterilisation systems:

① Lid	JM020
LxWxH 539 x 152 x 28 mm	
② Body	JM021
LxWxH 540 x 144 x 65 mm	
③ Single-use filter	JF167
LxW 235 x 118 mm	
④ Universal filter holder	JK098
LxW 245 x 100 x 6 mm	
⑤ Tray	JF435R
LxWxH 501 x 134 x 45 mm for up to two laparoscopic scopes and one fibre-optic cable	
⑥ Eye-piece holder	JG367
⑦ Shaft holder	JG368
⑧ Identification label	JG750
LxH 67 x 35 mm with adhesive back and process indicator	





STERRAD® is a registered trademark of Advanced Sterilization Products, a Johnson& Johnson company.
For more information on STERRAD®, refer to: www.aspjj.com

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
Phone (0 74 61) 95-0 | Fax (0 74 61) 95-26 00 | www.aesculap.de

Aesculap – a B. Braun company

Subject to technical alterations. This brochure is to be used exclusively for the offering, purchase and sale of our products. Reprinting, in whole or in part, prohibited. We reserve the right to withdraw the catalogues and price lists and to seek redress in the event of misuse.

Brochure No. C62802

1109/1.0/1